



Laser based technique for CO₂ flux measurements

Danijela Smajgl and Magda Mandic

Thermo Fisher Scientific, Bremen, Germany (danijela.smajgl@thermofisher.com)

The isotopic composition of CO₂ has been widely used for studying ecosystem gas exchange. The carbon isotope ratio and oxygen isotope ratio of atmospheric CO₂ ($\delta^{13}\text{C}$; $\delta^{18}\text{O}$) can be used to partition the gross fluxes of CO₂ in terrestrial ecosystems, such as plant respiration, soil respiration, and plant assimilation, as well as a tool for studying ocean-atmosphere interactions.

The characteristic $\delta^{13}\text{C}$ value is modified by plant metabolism and photosynthesis, while the $\delta^{18}\text{O}$ is affected by the oxygen exchange between the molecules of CO₂ and H₂O which originate from different water pools.

Here we present a new possibility for monitoring of CO₂ fluxes, Thermo ScientificTM Delta RayTM spectrometers. Instrumentation is based on direct absorption spectroscopy and it uses mid-infrared laser that operates at 4.3 μm . The laser scans over the spectral region containing four CO₂ absorption lines and isotope ratios are calculated from the spectrum fit. Calculation of different CO₂ isotopologues and determination of stable isotope ratios from spectrum is possible due to absorption lines which are shifted relative to each other. In direct absorption laser spectroscopy, the calculated isotope ratios are highly dependent on the CO₂ concentration. Therefore, for maximal accuracy of the measurements, instrument has a possibility to adjust the concentration of reference gas to match the concentration of the sample gas. This technology enables the simultaneous determination of $\delta^{13}\text{C}$, $\delta^{18}\text{O}$ and CO₂ concentration with a precision better than 0.05‰. With measurement resolution of seconds, it is a great tool for monitoring quick changes in open systems and provides new insights into CO₂ fluxes.