

Table-top Cathodoluminescence Microscopy for Geology



Toon Coenen

Business unit owner CL solutions coenen@delmic.com



Prof. Albert Polman **Dion Ursem** Hans Zeijlemaker

Thomas van der Heijden Sudiksh Srivastava Dorus elstgeest Wouter Roelofsen Andries Effting Sander den Hoedt

6th May 2020

Table of contents



- Slides 3 5: Floor model SEM CL
- Slides 6 8: Table-top CL
- Slide 9: More CL
- **Slide 10 16:** Supplementary slides

Cathodoluminescence system on floor model SEM





Floor model SEM-CL data examples

Ŧ

Floor-model SEM CL

Question: How do we make CL imaging more accesible and simpler?

+ High data quality

- High price
- + High spatial resolution
- + Modular and flexible
- + Correlation with other advanced SEM techniques
- Dedicated lab space and complex infrastructure needed
- Large footprint
- Expert user

Table-top CL

Table-top scanning electron microscopes are compact, user friendly, and affordable. Currently, there are no proper CL solutions for such systems, however. Here, we integrate a CL system on a table-top SEM

Cathodoluminescence

- + Low price
- + Simple to use
- + Fast
- + Small footprint
- + Little infrastructure required

- Lower spatial resolution
- Smaller range of acceleration voltages

System approaches

Fiber based collection system

Collection fiber & GRIN lens

Spectrometer

- Simple collection system with fiber and GRIN lens
- Direct coupling into spectrometer
- Relatively low collection
 efficiency due to reduced NA

- Higher collection efficiency
- Design is finished
- Results coming soon!

First results with fiber based CL system

- Hyperspectral CL maps for zircons with a full spectrum in every scanning pixel
- ➢ BSE and CL image show inverted contrast which is common in zircons
- Emission from Lanthanide ions is visible in the spectrum

More results are coming, stay tuned!

8

More CL

If you have any questions about cathodoluminescence imaging in general or about table-top cathodoluminescence feel free to contact me at <u>coenen@delmic.com</u> or visit our website www.delmic.com

More CL imaging at the EGU: Please visit this presentation if you are interested!

EGU2020-20478 Session BG4.4 "Correlative cathodoluminescence and EDS imaging of the benthic agglutinated foraminifer Liebusella goesi" Sangeetha Hari *et al.*, Wednesday, 06 May 2020, 08:30-10:15

Supplementary slides

Cathodoluminescence generation

Cathodoluminescence is the process whereby light (UV-VIS-IR) is generated when an electron beam hits a specimen. The emitted light carries a signature of the electronic structure of the material

Excited State

Ground State

Cathodoluminescence process in rocks

For a crystalline material, electrons in that material can only occupy certain energy states. Typically, (almost) all electrons reside in the valence band

- Rocks are typically insulators with wide band gaps between 5 15 eV (DUV-EUV)
- In CL we measure in the 0.8 6 eV range
- Defect states play an important role

Defect emission in rocks

CL versus other SEM-based techniques

CL can be used to extract various types of information and gives a unique contrast

Sample preparation for SEM-CL

CL imaging modes

- Measure CL intensity
- Short dwell times (10 100 µs) → videorate imaging
- Coarse spectral filtering and RGB mapping

- Measure CL spectrum
- Longer dwell times (10 1000 ms)
- Hyperspectral imaging with high spectral resolution