



## Automated CO<sub>2</sub> extraction from air for clumped isotope analysis in the atmo- and biosphere

Magdalena Hofmann (1), Martin Ziegler (1,2), Thijs Pons (1), Lucas Lourens (1), and Thomas Röckmann (1)  
(1) University of Utrecht, Utrecht, Netherlands, (2) ETH, Zurich, Switzerland

The conventional stable isotope ratios <sup>13</sup>C/<sup>12</sup>C and <sup>18</sup>O/<sup>16</sup>O in atmospheric CO<sub>2</sub> are a powerful tool for unraveling the global carbon cycle. In recent years, it has been suggested that the abundance of the very rare isotopologue <sup>13</sup>C<sup>18</sup>O<sup>16</sup>O on m/z 47 might be a promising tracer to complement conventional stable isotope analysis of atmospheric CO<sub>2</sub> [Affek and Eiler, 2006; Affek et al. 2007; Eiler and Schauble, 2004; Yeung et al., 2009]. Here we present an automated analytical system that is designed for clumped isotope analysis of atmo- and biospheric CO<sub>2</sub>.

The carbon dioxide gas is quantitatively extracted from about 1.5L of air (ATP). The automated stainless steel extraction and purification line consists of three main components: (i) a drying unit (a magnesium perchlorate unit and a cryogenic water trap), (ii) two CO<sub>2</sub> traps cooled with liquid nitrogen [Werner et al., 2001] and (iii) a GC column packed with Porapak Q that can be cooled with liquid nitrogen to -30°C during purification and heated up to 230°C in-between two extraction runs. After CO<sub>2</sub> extraction and purification, the CO<sub>2</sub> is automatically transferred to the mass spectrometer. Mass spectrometric analysis of the <sup>13</sup>C<sup>18</sup>O<sup>16</sup>O abundance is carried out in dual inlet mode on a MAT 253 mass spectrometer. Each analysis generally consists of 80 change-over-cycles. Three additional Faraday cups were added to the mass spectrometer for simultaneous analysis of the mass-to-charge ratios 44, 45, 46, 47, 48 and 49. The reproducibility for δ<sup>13</sup>C, δ<sup>18</sup>O and Δ<sub>47</sub> for repeated CO<sub>2</sub> extractions from air is in the range of 0.11‰ (SD), 0.18‰ (SD) and 0.02 (SD)‰ respectively.

This automated CO<sub>2</sub> extraction and purification system will be used to analyse the clumped isotopic signature in atmospheric CO<sub>2</sub> (tall tower, Cabauw, Netherlands) and to study the clumped isotopic fractionation during photosynthesis (leaf chamber experiments) and soil respiration.

### References

- Affek, H. P., Xu, X. & Eiler, J. M., *Geochim. Cosmochim. Acta* **71**, 5033–5043 (2007).  
Affek, H. P. & Eiler, J. M., *Geochim. Cosmochim. Acta* **70**, 1–12 (2006).  
Eiler, J. M. & Schauble, E., *Geochim. Cosmochim. Acta* **68**, 4767–4777 (2004).  
Werner, R.A., Rothe, M. & Brand, W.A., *Rapid Comm. Mass Spec.* **15(22)**, 2152-2167 (2001)  
Yeung, L. Y. *et al.*, *Proc. Natl. Acad. Sci.* **106**, 11496–11501 (2009).