



Phosphorus L_{2,3}-edge NEXAFS spectroscopy for the quantification of organic and inorganic P compounds in agro-industrial byproducts and soils

Luis Carlos Colucho Hurtarte (1), Helen Carla Amorim (2), Jens Kruse (3), Liming Wang (4), Julio Criginski Cezar (5), and Jörg Prietzel (1)

(1) Technische Universität München, Lehrstuhl für Bodenkunde, Freising, Germany (luis.colucho@wzw.tum.de), (2) Universidade Federal de Lavras, Campus Universitário, Lavras, MG 37200-000, Brazil, (3) Institute of Bio- and Geosciences, Agrosphere Institute (IBG-3), Forschungszentrum Jülich GmbH, Wilhelm Johnen Straße, 52425 Jülich, Germany, (4) Central Institute for Engineering, Electronics and Analytics, Analytics (ZEA-3), Jülich, Germany, (5) Laboratório Nacional de Luz Síncrotron – Centro Nacional de Pesquisa em Energia e Materiais – CP 6192, 13083-970 Campinas, SP, Brazil

Phosphorus (P) is an essential element for life on Earth that moves through its global cycle in different forms and associations. In contrast to inorganic P forms, organic P forms in soils and sediments have been historically unaccounted for due to the difficulty to accurately determine organic P forms in association with complex materials like soils and sediments. ³¹P NMR spectroscopy currently is the most reliable technique for bulk P speciation in soils and sediments. However, it does not account for all of the P present in soils and sediments, because alkaline extraction is required as pre-treatment, which does not completely recover soil- or sediment-bound P. During recent decades, synchrotron-based X-ray absorption spectroscopy (XAS) techniques (XANES, NEXAFS) have emerged as powerful tools for a direct, non-invasive P speciation in soils and sediments, permitting bulk and microspatial speciation analyses. Until now, most P XAS studies on soils and sediments have been conducted in the energy range of the P K-edge, which is well suited for the identification and quantification of different inorganic P species, but not for organic P species. XAS studies conducted at the P L-edge, which is a feature-rich edge for both organic and inorganic P forms are still scarce, because spectrum deconvolution is difficult. Thus, in this work we aim to apply XAS at the P L-edge for the quantification of different inorganic and organic P forms in soils and other environmentally important samples. First, we developed a data workflow to yield Linear Combination Fitting (LCF) results using known standard mixtures, then we analysed different soil samples and compared the P speciation obtained by XAS with results obtained by ³¹P NMR spectroscopy. Background subtraction and normalization were achieved by subtracting a polynomial background from the spectra and normalizing to the highest peak, using our R package “LCF”. The results obtained for standard mixtures with known composition showed that NEXAFS spectroscopy performed at the P L-edge has a high accuracy for the quantification of inorganic P (apatite), P-monoesters, P-diesters and phosphonates. However, it underestimates polyphosphates, especially when mixed with P-monoesters. The underestimation is possibly due to the large similarity between P-monoester and polyphosphate spectra. We also observed spectral distortions caused by the Si L-edge, which we could account for by subtracting an appropriate polynomial background to our data. To prove the feasibility of our technique for the P speciation of environmental samples, we analysed soil samples and agro-industry by-products and their EDTA-extracts by P L-edge NEXAFS, and compared the P speciation results with those obtained by ³¹P NMR spectroscopy. In conclusion, P L_{2,3}-edge NEXAFS spectroscopy is a potentially powerful tool for P speciation in geosciences, when combined with microscopy or as standalone spectroscopy.