

A new furnace design to improve the thermal demagnetization in paleomagnetism

Huafeng Qin, Xiang Zhao, Shuangchi Liu, Greig A. Paterson, Zhaoxia Jiang, Shuhui Cai, Qingsong Liu, Rixiang Zhu 2020

Content

- Introduction
- New design of furnace for thermal demagnetization
- Thermal demagnetization experiments
- Conclusions



The samples should be in zero magnetic field during thermal demagnetization. However, magnetic field noise, including residual magnetic fields of material and induced fields caused by the heating current in the furnace are always present.

Configuration of non-inductive heating wires





solenoids arranged in opposite direction (a), a bifilar solenoid wrapped (b) and new heating element made by solenoid and connected with a straight wire in series at one end for new developed demagnetizing furnace (c).

Simulation of solenoid magnetic field



Magnetic field measurement



Degaussing the shield cylinder



A solenoidal (a) and toroidal (b) degaussing coil embody the shield cylinder, and residual field in in shielding cylinder after demagnetization (c) with winding method according to figure (b). After demagnetization with this procedure, the residual magnetic field in the sample region can be less than 1 nT.

Temperature control



The heating current can be adjusted by Single-phase Voltage Module, which can reduce the current when approaching the target temperature. Its advantage is to ensure temperature stability and reduce the induced magnetic field. The temperature of heating process can be recorded in USB memory.

Experimental Method



Specimen and Demagnetization of TRM



Re-demagnetization of Synthetic specimens



It proved that the extra undesirable remanence acquired during the thermal re-demagnetization process

Reason for obtaining remanence



Nosie Remanece Vs NRM ('TRM', 'DRM')



Suppose that detrital remanence (DRM) is typically an order of magnitude weaker than TRM, nosie fields in furnace may have a large impact on thermal demagnetization data obtained from sedimentray rocks

Thermal demagnetization of natural specimens



Conclusions

- 1. The new heating element design can significantly reduce magnetic field noises in thermal demagnetizer
- 2. The experiments confirm that magnetic field noise in various furnaces can have an observable and detrimental impact on demagnetization behavior
- 3. The new furnace represents an important improvement in the design of thermal demagnetizers and allows for extremely weak specimens to be successfully measured.

Acknowledgements

We thank Jie Yuan and Dr Shihu Li for providing the natural samples for testing. This work was supported by NSFC grants 41674073 and the Project of National Deep Exploration Technology and Experimental Research(Grant number:SinoProbe-09-02(201011079)).

