

## ”ALPOL” experiment on the Russian Luna-Globe and Lunar-Resource missions

**M.V. Gerasimov** (1), I.I. Vinogradov (1), A.G. Sapgir (1), A.V. Kalyuzhnyi (1), O.E. Kozlov (1), A.V. Stepanov (1), A.Yu. Titov (1), D.S. Rodionov (1), R.S. Zverev (1), M.M. Gerasimov (1), M. Cabane (2), C. Szopa (2), D. Coscia (2), J.-P. Goutail (2) J.J. Correia (2), A. Galic (2), F. Goesmann (3), P. Wurz (4), D. Abplanalp (4), and M. Tulej (4),  
(1) Space Research Institute of the RAS, Moscow, Russia, (2) LATMOS, IPSL, CNRS, Paris, France (3) Max Planck Ins. for Solar System Res., Katlenburg-Lindau, Germany, (4) Physics Institute, University of Bern, Switzerland;  
(mgerasim@mx.iki.rssi.ru)

### Abstract

The paper considers possible sources of polar volatile depositions on the Moon and presents the experiment which is in preparation for Russian polar landers to characterize polar regolith volatiles.

### Introduction:

Luna-Globe and Lunar-Resource missions are now under preparation in Russia for launch in 2015 (16?). Both missions have landers to be put in northern and southern polar regions of the Moon for investigation of possible depositions of water and other volatile components. The gas-analytic experiment “ALPOL” is aimed on comprehensive investigation of the inventory of volatiles in the regolith of the lunar polar regions. This information is valuable for understanding of the real content and chemical composition of volatiles in polar regions regolith. Chemical composition and abundance of volatile components bears information about the source and dynamics of volatiles on the Moon. Information about isotopic ratios of key volatile elements (C, H, O, N, noble gases, ...) is also of great importance for understanding of sources of the lunar volatile components.

### Tasks of the Gas-Analytic Package:

The main tasks of the Gas-Analytic Package (GAP) are:

1. Detailed investigation of chemical composition and abundances of volatile compounds (H<sub>2</sub>O, CO<sub>2</sub>, N<sub>2</sub>, H<sub>2</sub>, noble gases, organics, etc.) in the surface regolith material of the Moon at the landing place;
2. Investigation of forms of incorporation of volatile components into the solid surface

material;

3. Investigation of organic components in the surface material;
4. Measurement of isotopic composition of CHON elements (<sup>13</sup>C/<sup>12</sup>C, D/H, <sup>17</sup>O/<sup>16</sup>O, <sup>18</sup>O/<sup>16</sup>O, <sup>15</sup>N/<sup>14</sup>N) and noble gases.
5. Measurements of the exosphere composition of the Moon.

### Description of the GAP:

Design of the GAP is based mainly on the heritage of the GAP developed for the Russian “Phobos-Grunt” mission [1, 3]. The GAP consists of three individual instruments: 1) Thermal Analyzer (TA-L); 2) Gas Chromatograph (GC-L); and 3) Neutral Gas Mass-Spectrometer (NGMS).

### TA-L instrument

TA-L instrument is aimed to release volatiles from the lunar regolith and to measure temperatures of the release. The volatiles can be frozen on regolith particles surfaces and/or bound to particles matrixes. The release of frozen volatile components can occur at a very low temperature range while the bound volatiles are released at rather high temperatures. For this reason TA-L has two types of ovens. The first one is a multi-use low temperature (<250°C) oven which is sealed immediately after direct load from the manipulator. Other ovens are a single-use high-temperature ones (~1000°C). TA-L has the Soil Acquisition System (SAS) and a carousel with 8 single-use pyrolytic cells (PC). SAS receives a portion of regolith from the manipulator and loads the high-temperature PC. PCs provide programmed heating of the imbedded regolith sample. Released gases are transferred to GC-L instrument via transfer capillary tubes.

Tasks of the TA-L instrument are:

1. To release and transfer to GC-L frozen and bound volatile components.
2. To measure exo- and endothermal reactions in the soil sample to determine minerals with phase transitions at temperatures < 1000°C;
3. To perform pyrolysis of heavy organics (kerogens?) and provide their transfer into GC-L and NGMS instruments.

## GC-L instrument

GC-L has two capillary columns (CC): one coated by mole sieves 5A for analysis of permanent gases; and the second coated by Porapac-Q for analysis of high boiling components. Each column is equipped by thermal conductivity detector (TCD). GC-L instrument has two injection traps (IT), the first filled by Carbosieve and the second by Porapac-Q absorbents, which can be cooled by Peltier elements down to -50°C. ITs sequentially collect high boiling components and permanent gases for injection by pulse heating into respective CCs.

GC-L instrument is a combination of a gas chromatograph with a tunable diode laser absorption spectrometer (TDLAS). TDLAS has four tunable diode lasers aimed on the measurement of CO<sub>2</sub> and H<sub>2</sub>O molecules and the C, O, and H isotopes in these molecules.

Tasks of the GC-L instrument are:

1. Accumulation of gases that are released from the regolith sample during pyrolysis;
2. Redistribution of gases of different types (permanent gases, organics, etc.) between respective columns;
3. Measurement of the abundance of each separate gas component;
4. Measurement of the isotopic ratios of <sup>13</sup>C/<sup>12</sup>C, D/H, <sup>17</sup>O/<sup>16</sup>O, <sup>18</sup>O/<sup>16</sup>O, in CO<sub>2</sub> and H<sub>2</sub>O.

## NGMS instrument

The NGMS instrument is a time-of-flight mass-spectrometer based on an earlier development for stratospheric research [2].

Tasks of the NGMS instrument are:

1. Mass spectrometric identification of gas components that are released from the gas chromatograph;
2. Measurement of elements isotopic ratios.
3. Measurements of the exosphere composition.

## Method of analysis.

The GAP receives a portion of soil from the Sampling Device of the manipulator. This portion is loaded by decision in one case directly into low-temperature oven or in the second case into the SAS of the TA-L instrument. Low-temperature oven provides release of frozen components and TA-L transfers them into GC-L. The SAS performs selection of a dose of the regolith sample (~0.5 cm<sup>3</sup>), loads and seals it in the high-temperature PC. The PC performs programmed heating of the sample up to 1000°C to do thermal analysis of the sample and provides the release of volatiles into gaseous phase. Released gases are transported to the gas chromatograph GC-L by a flow of carrier gas (helium). Gases in time of pyrolysis are analyzed in GC-L using TDLAS to measure H<sub>2</sub>O and CO<sub>2</sub> molecules and isotopic ratios <sup>13</sup>C/<sup>12</sup>C, D/H, <sup>17</sup>O/<sup>16</sup>O, <sup>18</sup>O/<sup>16</sup>O in them. Afterwards gases are separated between two injection traps and accumulated in them. Collected gases are then analyzed on both chromatographic capillary columns. Separated gases are transferred to mass spectrometer for mass spectrometric analysis for their proper identification and measurement of isotopic composition of key volatile element.

## Cooperation:

Main partners of the Gas Analytic Package team are: Space Research Institute (IKI) of the RAS (Russia), Vernadsky Inst. of Geochem. (GEOKHI) (Russia), LATMOS and LISA, University of Paris (France), Max Planck Ins. for Solar System Res. (Germany) Physics Institute, University of Bern, (Switzerland)

## References

- [1] M. Gerasimov, et al. Gas-chromatographic complex (TDA, KHMS-1F, MAL-1F), In: Phobos-Grunt, the project of space expedition, ed.: A.V. Zakharov, ed. Lavochkin Association and IKI, pp. 249-270, 2011.
- [2] D. Abplanalp, P. Wurz, L. Huber, I. Leya, E. Kopp, U. Rohner, M. Wieser, L. Kalla, and S. Barabash, "A neutral gas mass spectrometer to measure the chemical composition of the stratosphere," *Adv. Space Res.* 44, 870–878, 2009.
- [3] L.M. Zelenyi, A.V. Zakharov, G.M. Polischuk, and M.B. Martynov, Project of the Mission to Phobos, *Solar System Research*, Vol. 44, 15–25, 2010.